

FITTING INSTRUCTIONS



FOR THE ELBOW EXTENSION

1. Place the patient in a reclined or resting position.
2. ONLY open FRONT CUFFS #1 & #2.
3. Position patient's arm with thumb in neutral or thumbs up position.
4. Line cams up with the elbow joint.
5. Close upper cuff #1 (FRONT OF UPPERARM) and pull to line.
6. Bend splint up at hinge to meet bottom of forearm.
7. Close cuff #2 (TOP OF FOREARM) and pull to line.
8. Support the arm with pillows to elevate the arm.

OPTIONAL: If Resting Hand Wrist Orthosis is being worn in conjunction with the Elbow Extension Dynasplint® System, apply the orthosis first, following the enclosed instructions in packaging.

Stretch Beyond Your Expectations.®

Taller patients will be fit with two - 2" straps in place of butterfly strap



DO NOT OPEN CUFFS THAT HAVE BEEN SECURED WITH DUCT TAPE OR GLUE