

## SERVICING THE CONTRACTURE MANAGEMENT NEEDS OF SKILLED NURSING FACILITIES

CONTACT US — Email: [helpingsnfs@dynasplint.com](mailto:helpingsnfs@dynasplint.com) | Fax: 866.887.2050 | Phone: 866.825.6179

### NeuroFlex™ Restorative™ Flex Hand (20721)

L Code: L3807

- Features:
- Wrist and palmar bar flex to correct wrist flexion.
  - Progressive sized hand cones or cylinders provide extension of the MP and IP joints.
  - Includes: Regular Cones/Optional Reverse Cones/Optional Cylinders



### BendEase Hand (20715)

L Code: L3807

- Features:
- Foam over aluminum technology is used to comfortably & quickly address lost ROM.
  - Easily shaped to conform to the patient's wrist/hand.

Size	Width of MCP	Wrist to Finger	
Small - Right or Left	Up to 3"	Up to 6 1/2"	<i>Universal width of palm at MCPs Length between wrist and fingertips right or left extremity</i>
Medium - Right or Left	Up to 3 1/2"	Up to 7"	
Large - Right or Left	Up to 3 3/4"	Up to 8"	



### Progressive Hand DHXY

L Code: L3924

- Features:
- Foam over aluminum technology is used to comfortably & quickly address lost ROM.
  - Custom fit specific to your patients' needs.

Use sizing chart  
for BendEase Hand, above.



### NeuroFlex™ Restorative™ Elbow (20724)

L Code: L3760

- Features:
- Universal extremity, hinges adjust to gradually increase extension of joint.
  - Lower cuff swivels to allow for changing angles of the extremity.
  - Double padded for skin protection.

	<b>Top:</b> Measure at lateral aspect - 1/2 of distance from shoulder joint to elbow joint*	<b>Bottom:</b> Elbow joint to approximately 1/2 the distance to hand MCP joints	Upper Arm Circumference	Lower Arm Circumference	
	3 1/2"	4"	6 1/2"-8 1/2"	5"-7"	<i>Child</i>
	4"	5"	8"-10"	6 1/2"-8 1/2"	<i>Ex-Small</i>
*Allow 2 finger width distance below axilla.	5"	6"	8"-12"	7"-10"	<i>Small</i>
	6"	6"	8"-13"	7 1/2"-11"	<i>Medium</i>
	6 1/2"	6 1/2"	10"-15"	8"-12"	<i>Large</i>



### RestAir™ Hip Orthosis (20760air)

L Code: L1652

- Features:
- Will address severe hip adduction and/or knee scissoring.
  - Abduct hip safely and comfortably.
  - Allows air flow for hygiene.
  - Inflate with same number of pumps in each of the 4 air bladders.
  - Spreader Bar Assembly can be added.

Thigh Circumference	
6"-12"	<i>Ex-Small</i>
11"-17"	<i>Small</i>
16"-22"	<i>Large</i>
21"-26"	<i>Ex-Large</i>



### NeuroFlex™ Restorative™ Knee (20725)

L Code: L1831

- Features:
- Universal extremity, hinge adjusts to gradually increase extension of joint.
  - Lower cuff swivels to allow for changing angles as the joint is corrected.
  - Double padded for skin protection.

Knee joint to approx mid femur	Knee joint to approx 1/2 the distance to lateral malleolus	Thigh Circumference	Calf Circumference	
5"	6"	8"-12"	7"-10"	<i>Child</i>
6"	6"	8"-13"	7 1/2"-11"	<i>Ex-Small</i>
6 1/2"	6 1/2"	10"-15"	8"-12"	<i>Small</i>
7 1/2"	7 1/2"	13"-17"	10"-15"	<i>Medium</i>
8 1/2"	8 1/2"	16"-22"	14"-18"	<i>Large</i>



### RM Safe Boot II (20744)

L Code: L4396/L4397

- Features:
- Non-skid sole attaches with hook and loop.
  - Design helps to maintain foot alignment while flex addresses the tone.

Foot Length	Youth Shoe	Men's Shoe	Women's Shoe	
Up to 8"	1.5		4	<i>Small</i>
Up to 9 1/2"		6.5	8	<i>Medium</i>
Up to 11"		11	12	<i>Large</i>
Up to 13"		16		<i>Ex-Large</i>

